

KeepKleen® 5 mil.

# DISPOSABLE LATEX

## Medical Exam Grade Gloves

### FUNCTIONAL BENEFITS:

Gloves are Class-2 medical grade and CFIA-approved for food handling

Manufactured in an ISO 9002 registered facility

5 mil. thickness and rolled cuffs offers greater tear resistance

Choice of white or blue. White available as powder-free

Excellent touch-sensitivity and comfort

### APPLICATIONS:

Dental and Medical Exam  
Food Processing  
Laboratories  
Pharmaceutical  
Janitorial



RDL

RDLPF

### KeepKleen® 5 mil. Disposable Latex Medical-exam Grade Gloves

Superior's KeepKleen® Class-2, medical-grade, disposable latex gloves are manufactured in an ISO 9002 registered facility. Recognized by the Canadian General Standards Board, they are certified and tested under ISO 11193.

At a thickness of 5 mil., they offer greater resistance to tearing, while providing superior touch-sensitivity, dexterity and comfort. Gloves have a beaded cuff-edge for extra strength that prevents ripping during donning.

Two colors are available: blue and white. White and blue versions are lightly powdered with food-grade cornstarch to ease donning and doffing. White is also available as a powder-free glove. All are CFIA-approved for food handling.

Caution: Latex gloves may cause allergic reactions in some individuals.

For more information on Disposable Gloves click the following link:

[http://www.superiorglove.com/Disposable\\_Rubber\\_Gloves\\_S26.html](http://www.superiorglove.com/Disposable_Rubber_Gloves_S26.html)

- |        |  |
|--------|--|
| RDL    | KeepKleen® 5 mil, White latex, lightly powdered, medical-exam grade, textured grip.<br>CFIA-approved for food handling |
| RDLPF  | KeepKleen® 5 mil, Powder-free, white latex, medical-exam grade, textured grip.<br>CFIA-approved for food handling      |
| RDLSPF | KeepKleen® 5 mil. latex, powder-free, smooth surface<br>Sizes: XS, S, M, L, XL   |
| RDLB   | KeepKleen® blue, 5 mil. latex, lightly powdered<br>Sizes: S, M, L, XL  |



LATEX Allergy



Use a NON-LATEX glove.

### Symptoms of Latex Allergies

- Red, crusty rash stopping at the wrist where the glove ends
- Red, itchy bumps sometimes associated with runny nose, red eyes and/or asthma
- Rare and severe latex allergies can even result in shock

### How to Avoid Latex Allergies

- Latex allergies are often aggravated by powdered gloves which can carry latex proteins. Use a non-powdered glove.
- Best of all, use a non-latex glove such as a nitrile or vinyl.

KeepKleen



Engineered hand protection for every industry®

### For Ordering Info or Questions:

Phone (519) 853-1920

Fax (519) 853-4496

Toll Free (800) 265-7617

Email [sales@superiorglove.com](mailto:sales@superiorglove.com)

[www.superiorglove.com](http://www.superiorglove.com)